

FOXTR 1



**FOXTR**  
ULTRA SUSPENSION



## FOXTR 1 Travel Mobility Scooter



The **FOXTR 1 Travel Mobility Scooter** with Ultra Suspension offers all the features you need to regain your independence in both indoors and outdoors spaces. This feature packed travel scooter offers larger wheels for a smoother ride, class-leading shock absorption, quick disassembly and a powerful motor.

Weight Capacity

**300 lbs**  
133 kg

Max. Range

**10 miles**  
16 km

Max. Speed

**4.4 mph**  
7 km/h

Turning Radius

**51"**  
129 cm

# FOXTR 1 Travel Scooter

**Weight Capacity:** 300 lbs (136 kg)

**Maximum Range:** 10 miles (16 km)

**Maximum Speed:** 4.4 mph (7 km/h)

**Climbing Angle:** 8 degrees

**Ground Clearance:** 3.1" (7.9 cm)

**Turning Radius:** 51" (129 cm)

**Dimensions:** 44.5" (L) x 22" (W)

**Seating:** 18" (W) x 16" (D) Padded, Sliding/Pivoting Seat

**Front Tires:** 2.1" x 7.5" Solid

**Rear Tires:** 2.8" x 8.5" Solid

**Batteries:** 12AH x12V (2)

**Heaviest Piece:** 39 lbs (17.7 kg)

**Battery Weight:** 20 lbs pair (9 kg)



## Ultra Suspension System

Advanced Front and Rear Suspension for improved outdoors performance



## Larger Solid Tires

Treaded 7.5" front and 8.5" back tires are larger than other travel scooters

## Features & Benefits

- Larger wheels than other travel scooters for a smoother ride.
- Forward and backward adjustable for ideal positioning
- Larger shopping basket
- Powerful, front LED headlight
- Ultra Suspension with front and rear shock absorbers
- Front and rear mudguards
- Disassembles into 5 separate parts for quick transportation
- Swiveling, adjustable-height seat
- Flip-up, width-adjustable armrests
- Tiller can be easily adjusted with a comfortable grip knob
- Durable rubber floor mat and foot rests
- Includes anti-tipping wheels & double mirrors for safety

**Warranty:** Limited lifetime warranty on the frame, 24 months on Electronics and 12 months on the Batteries



## Intuitive Control Panel

Easy-read battery gauge, speed control, horn and lights



## Adjustable Seating

Folding seat can pivot and slide to help you find perfect fit



## Removable Battery Pack

Removable batteries can be charged onboard scooter or separately