

FOXTR 3



FOXTR3  
ULTRA SUSPENSION



## FOXTR 3 Heavy Duty Mobility Scooter



The **FOXTR 3 Heavy Duty 4-Wheel Scooter** features sturdy, heavy duty construction to provide an exhilarating & comfortable ride. With an impressive range, advanced suspension system and an impressive suite of comfort features, the FOXTR 3 is ideal for all your outdoor adventures.

Weight Capacity

**425 lbs**  
193 kg

Max. Range

**25 miles**  
40 km

Max. Speed

**8 mph**  
13 km/h

Turning Radius

**53"**  
134 cm

# FOXTR 3 Heavy Duty Scooter

**Weight Capacity:** 425 lbs (193 kg)

**Maximum Range:** 25 miles (40 km)

**Maximum Speed:** 8 mph (13 km/h)

**Climbing Angle:** 8 degrees

**Ground Clearance:** 4" (10.1 cm)

**Turning Radius:** 53" (134 cm)

**Dimensions:** 56" (L) x 26" (W) x 53" (H)

**Seating:** 20" (W) x 20" (D) x 22" (H) Captain's Chair

**Front Tires:** 14" x 4" pneumatic

**Rear Tires:** 14" x 4" pneumatic

**Batteries:** 50AH x12V (2)

**Weight without Batteries:** 185 lbs (84 kg)

**Battery Weight:** 32 lbs each (14.5 kg)



## Ultra Suspension System

Advanced Front and Rear Suspension for excellent outdoors performance

## Features & Benefits

- Full suspension and large, 14" pneumatic tires (air filled)
- Hazard lights, LED lighting package, turn signals & reflectors
- One-touch, infinitely adjustable delta tiller for comfort
- High-back captain's seat folds-down and is height adjustable
- Seat features reclining backrest and adjustable headrest
- Dual rear-view mirrors, anti-tip wheels & positioning belt
- Convenient USB port for charging your mobile device
- Large front plastic basket to carry personal belongings
- Wider 20" seat for comfort and versatility
- Includes dual mirrors for improved visibility

**Warranty:** Limited lifetime warranty on the frame, 24 months on Electronics and 12 months on the Batteries



## Large Pneumatic Tires

14" air-filled tires help you navigate challenging terrains in comfort



## Includes Dual Mirrors

Left and right mirrors included for better visibility



## Luxury Seating

Wider 20" Captain's Chair reclines, pivots & slides back & forth



## Premium Lighting Package

Hazard lights, LED lighting, turn signals and reflectors