

FOXTR 2



FOXTR2
ULTRA SUSPENSION



FOXTR 2 Mid-Size Mobility Scooter



The FOXTR 2 Mid-Size Mobility Scooter incorporates innovative design to guarantee a comfortable ride. With a luxurious captain's seat and large solid tires, the FOXTR 2 excels at navigating diverse terrains.

Weight Capacity

300 lbs
181 kg

Max. Range

20 miles
32 km

Max. Speed

5.7 mph
9.2 km/h

Turning Radius

60"
152 cm

FOXTR 2 Mid-Size Scooter

Weight Capacity: 300 lbs (181 kg)

Maximum Range: 20 miles (32 km)

Maximum Speed: 5.7 mph (9.5 km/h)

Climbing Angle: 12 degrees

Ground Clearance: 3.5" (8.9 cm)

Turning Radius: 60" (152 cm)

Dimensions: 51" L x 25" W

Seating: 19"W x 16"D x 21"H - Captain's Chair


Front Tires: 3"x 10" solid

Rear Tires: 3"x 10" solid

Batteries: 38AH x12V (2)

Weight without Batteries: 194 lbs (88 kg)

Battery Weight: 24 lbs each (10.8 kg)



Ultra Suspension System

Front and Rear
Suspension for a more
comfortable indoors /
outdoors ride



Full Lighting Package

Includes hazard lights, LED
lighting package, turn
signals and reflectors

Features & Benefits

- Ultra suspension and large, 10" solid tires
- Includes hazard lights & LED lighting package
- Turn signals and reflectors for added safety
- One-touch, infinitely adjustable delta tiller
- High-back captain's seat folds-down & is height adjustable
- Includes dual rear-view mirrors
- Anti-tip wheels, free-wheel mode & positioning belt
- Large front plastic basket to carry your personal belongings
- Two Colour Options: Available in Chili Red and Sky Blue

Warranty: Limited lifetime warranty on the frame, 24 months
on Electronics and 12 months on the Batteries



Includes Dual Mirrors

Left and right
mirrors included
for better visibility



Large Solid Tires

10" solid tires for
reliability and comfort



Roomy Carry Basket

For everyday errands and storage
of personal belongings