

VITAL-WH18



Vital Lightweight Folding Wheelchair



The **Vital Lightweight Wheelchair** is a highly portable, sturdy manual wheelchair perfect for home use and errands around your community. It folds compactly to fit even the tightest storage spaces. It includes easy-flip back padded full arms, swing-away footrests, and hand brakes for convenient attendant care.

Weight Capacity

220 lbs
100 kg

Seat Width

18"
45.7 cm

Seat Depth

16"
40.6 cm

Footrests

**Swing
Away**

Vital Lightweight Wheelchair

Weight Capacity: 220 lbs (100 kg)

Seat Dimensions: 18" x 16"

Overall Weight: 32 lbs (14.7 kg)

Seat-to-floor height: 19" (48.3 cm)

Overall Dimensions Open: 36" (L) x 24.5" (W) x 35"(H)

Overall Dimensions Closed: 36" (L) x 12.7" (W)

Rear Wheel Size: 22"

Front Caster Size: 7"

Carton Shipping Weight: 38 lbs (17.2 kg)

Shipping Dimensions: 30.7" x 13.4" x 32.2"

Footrests: Swingaway

Armrests: Full length, flip-back



Ultra Foldable Frame

Wheelchair folds in half and backrest folds back for ultimate portability

Features & Benefits

- Rust-resistant, powder-coated aluminum alloy frame
- Self-locking nursing break and wheel locks for safety & convenience
- The wheelchair features a foldable backrest for optimal storage and transportation
- Flip-up full-length padded armrests
- Comfortable cushion crafted from oxford and breathable mesh fabric
- Detachable swing-away footrests
- Ergonomic hand rim is easy on your hands
- Double-cross axle for durability and support


Warranty: Limited lifetime warranty on frame.

90-day warranty on brakes, wheels, casters and other components.



Flip Back Armrests

Full, padded armrests for support can flip-back for lateral transfers and compatibility with tables and desks



Dual Brake System

Ergonomic hand brakes for caregiver convenience and easy to reach wheel locks



Quick-release footrests

Trigger release, length-adjustable, swing-away footrests



Ergonomic 22" Wheels

Molded hand rim for easier maneuverability and reduced fatigue