

VITAL-TC18



Vital Lightweight Folding Transport Chair



The **Vital Lightweight Transport Chair** is a compact and feature-packed transport wheelchair perfect for helping a loved one or client move around their home, enjoy outdoor time, and visit stores and doctors' offices. It folds compactly to fit in the tightest storage spaces. It includes easy-flip back padded full arms, a revolutionary flip back, adjustable footrests, and both foot and hand brakes for convenient attendant care.

Weight Capacity

220 lbs
100 kg

Seat Width

18"
45.7 cm

Seat Depth

16"
40.6 cm

Footrests

Flip Back

Vital Lightweight Transport Chair

Weight Capacity: 220 lbs (100 kg)

Seat Dimensions: 18" x 16"

Overall Weight: 25 lbs (11.3 kg)

Seat-to-floor height: 18.5" (46.9 cm)

Overall Dimensions Open: 36" (L) x 24.5" (W) x 35"(H)

Overall Dimensions Closed: 36" (L) x 12" (W)

Rear Caster Size: 16"

Front Caster Size: 7"

Carton Shipping Weight: 31 lbs (14 kg)

Shipping Dimensions: 28.7" x 11.8" x 31.5"

Footrests: Flip-Back

Armrests: Full-length, flip-back



Ultra Foldable Frame

Wheelchair folds in half and backrest folds back for ultimate portability

Features & Benefits

- Rust-resistant, powder-coated aluminum alloy frame
- Self-locking nursing break and wheel locks for safety & convenience
- Transport Chair features a foldable backrest for optimal storage and transportation
- Flip-up full-length padded armrests
- Comfortable cushion crafted from oxford and breathable mesh fabric
- Detachable, revolutionary flip back footrests
- Double-cross axle for durability and support

Warranty: Limited lifetime warranty on frame.

90-day warranty on brakes, wheels, casters and other components.



Flip Back Armrests

Full, padded armrests for support can flip-back for lateral transfers and compatibility with tables and desks



Dual Brake System

Ergonomic hand brakes for caregiver convenience and easy to reach wheel locks



Flip Back Footrests

Save space with the Vital Transport Chair's convenient flip back footrests



Larger 16" Rear Casters

Larger, treaded rear casters for excellent indoors/outdoors performance