

MultiTech



Technical specifications:

Width of available platforms	36" (89 cm)
	39" (96 cm)
	48" (120 cm)

Platform length:	75" or 80" (190 or 203 cm)
------------------	----------------------------

Minimum height:	12" (30,5 cm): Individual brake base
	12 ¼" (31 cm): Central brake base

Maximum height:	29 ½" (75 cm): Individual brake base
	29 ¾" (75.5 cm) : Central brake base

Load capacity:	500 lb (227 kg)
----------------	-----------------

Maximum backrest angle:	80°
-------------------------	-----

Maximum footrest angle:	45°
-------------------------	-----

Clearance under the bed:	5 5/8" (14 cm): Individual brake base
	5 ½" (18 cm): Central brake base

Bed weight:	209 lb (95 kg)
-------------	----------------

Features

- Bed designed for home care patients or long-term care facilities
- Ultra-quiet DC motors. Helps maintain a calm and relaxed environment
- Easy to maintain with quick connect components
- IP54-approved electrical components protected from dust and water infiltration
- Battery ensuring temporary autonomous operation
- Bed can be moved at any heights even when in its lowest position
- Antimicrobial and easy to clean powder coated steel structure
- 4 Casters equipped with individual brakes
- 4 " (10 cm) diameter medical casters
- 6 functions remote control
- Manually adjustable foot section angle (choice of two positions: flexion and standard)

Options

- Digital nurse control integrated in the footboard (include the head and footboard)
- Digital patient control integrated in the siderails (requires the nurse control, include one pair of siderail)
- Central brake system (not available on the 48" bed)

Accessories

- Composite (grey or beige) or melamine decorative headboard and footboard (Choice of melamine colors: white, rosewood, maple, and other colors on request)
- Composite platform cover
- Rotec tilt-up half side rails 28" (71 cm) (impossible to fit two pairs of this model)
- Rotec tilt-up half side rails 23" (58 cm)
- Rotec tilt-up half side rails 21" (53 cm) grey color
- Rotec full length side rails
- Metallic half side rails
- Rotec rotating half side rails 28" (71cm)
- Rotec rotating half side rails with increased protection height 28" (71cm)
- Rotec rotating half side rails 21" (71cm) grey color
- Bed length extension to 84" (213 cm)
- Extensions to increase the underbed clearance of 1" or of 2" (2,5 cm or 5 cm) for casters with individual brake
- Fixed position grab bar (unavailable on 48" bed)
- IV pole
- Bed trolley
- Metallic support with bar and trapeze
- Pump support
- Siderail pad
- Medical mattress

