



MultiTech ULB

Ultra-Low Bed



Technical specifications

Platform width: 35" (89 cm)

Platform length: 78" (198 cm)

Minimum height: 8" (21 cm)

Maximum height: 29" (74 cm)

Load capacity: 580 lb (264 kg)

Maximum backrest angle: 68° (degrés)

Maximum footrest angle: 34° (degrés)

Trendelenburg Position: -18° to +18° (degrés)

Dégagement sous le lit :

Variable according to the height of the bed:
6" (12 cm) clearance when bed is 17" high (43 cm)

Bed weight: 245 lb (111 kg)

Multitech ULB

Features

- Bed designed for home care patients or long-term care facilities
- Helps prevent Falls Related Injuries. Ultra-Low bed 8" (21 cm)
- Ultra-quiet DC motors. Helps maintain a calm and relaxed environment
- Battery ensuring temporary autonomous operation
- IP54-approved electrical components protected from dust and water infiltration
- Bed can be moved at any heights even when in its lowest position
- Antimicrobial and easy to clean powder coated steel structure
- Casters equipped with a synchronized braking system (2 x 2 casters)
- 4" (10 cm) diameter medical casters
- 8 functions remote control
- Manually adjustable foot section angle (choice of three positions: flexion, standard and extension)
- Side and foot mattress retainer
- 6 receptacles for IV pole
- Bumpers at the head and foot end of the bed
- Bed may be separated into two sections to facilitate transport
- Authentic Canadian product

Options

- Digital nurse control integrated in the footboard (include the head and footboard)
- Digital patient control integrated in the siderails (requires the nurse control, include one pair of siderail)

Accessories

- Composite or melamine decorative headboard and footboard (Choice of melamine colors: white, rosewood, maple, and other colors on request)
- Composite platform cover
- Rotec tilt-up half side rails 28" (71 cm)
- Rotec full side rails
- Rotec rotating half side rails 28" (71 cm)
- Rotec rotating half side rails with increased protection height 28" (71 cm)
- IV pole
- Bed trolley
- Bar with trapeze
- Siderail pad
- Pump support
- Medical mattress

