



RT-500



The RotaBed Rotating Chair Bed



The RotaBed Rotating Chair Bed is a game-changer in the world of home care and long-term care. This revolutionary rotating hospital bed allows individuals to get in and out of bed without any assistance at the simple touch of a button.



The RotaBed Rotating Chair Bed

Product Code: RT-500

Total Load capacity: 390 lbs / 175 kg

Min.Height: 15 inches / 38.5 cm

Max Height: 32 inches / 81 cm

Max. Backrest Angle: 70°

Max. Legrest Angle: 20°

Sleep Deck Dimensions: 35.4" x: 74.8" / 90 cm x 190 cm

Overall Dimensions: 40" x 84" / 102 cm x 213 cm

Trendelenburg Position: BOTH WAYS

Trapeze: OPTIONAL



Watch RotaBed Video



Efficient Mobility

Transitions between lay-down and seated positions in only 40 seconds



Safety First

Includes full locking safety sides and handrails for a safe patient experience

Features & Benefits

- Bed-to-Chair function rotates sleep surface to allow unassisted transfer to and from bed
- Goes from lying down to seated position in only 40 seconds
- Intelligent design automatically adjusts bed to correct height when bed-to-chair function is activated
- Easy to operate remote control offers endless positional configurations
- Includes Trendelenburg and Reverse Trendelenburg positions
- Optional Trapeze
- Equipped with locking safety sides and hand rails
- Beechwood headboard, footboard and the safety side-rails match most decor styles
- Requires articulating pressure prevention mattress

Warranty: 5 year warranty on the motor and limited lifetime warranty on the frame



Smart Technology

Bed height automatically adjusts for Bed-to-Chair function



Adaptable Comfort

Articulating mattress adapts to the RotaBed's movement



Ultimate Adjustability

Instinctive, ergonomic hand pendant offers endless positional configurations