

VariTech



Technical specifications

Platform width:	Adjustable from 35" to 38" to 47" to 53" (89 cm to 97 cm to 119 cm to 134 cm)	Maximum backrest angle:	64° (degrees)
Platform length:	78" (198 cm)	Maximum footrest angle:	36° (degrees)
Minimum height:	13" (33 cm) or 13 ½" (34 cm) for fixed base	Trendelenburg Position:	-14° to +14° (degrees)
Maximum height:	30" (76 cm) or 30 ½" (77 cm) for fixed base	Dégagement sous le lit:	Variable according to the height of the bed: 6" (12 cm) clearance when bed is 17" high (43 cm)
Load capacity:	1100 lb (499 kg)	Bed weight:	309 lb (140 kg)

Features

- Bed designed for home care patients or long-term care facilities
- Bed with the highest load capacity in its category
- Affordable width-adjustable bariatric bed
- Ultra quiet DC motors. Helps maintain a calm and relaxed environment.
- Width-adjustable bed platform that does not require tools
- Battery ensuring temporary autonomous operation
- IP54-approved electrical components protected from dust and water infiltration
- Bed can be moved at any heights even when in its lowest position
- Composite platform spring cover
- Antimicrobial and easy to clean powder coated steel structure
- Casters equipped with a synchronized braking system (2 x 2 casters)
- 4" / 10 cm diameter medical casters
- 8 functions remote control
- Manually adjustable foot section angle (choice of three positions: flexion, standard and extension)
- Side and foot mattress retainer
- 6 receptacles for IV pole
- Bumpers at the head and foot end of the bed
- Bed may be separated into two sections to facilitate transport
- Authentic Canadian product

Options

- Digital nurse control integrated in the footboard (include the head and footboard)
- Digital patient control integrated in the siderails (requires the nurse control, include one pair of siderail)
- Length extension of the platform to 88 ½" (225 cm)
- Fixed base

Accessories

- Composite or melamine decorative headboard and footboard (Choice of melamine colors: white, rosewood, maple, and other colors on request)
- Rotec tilt-up half side rails 28" (71 cm)
- Rotec rotating half side rails 28" (71 cm)
- Rotec rotating half side rails with increased protection height 28" (71 cm)
- IV pole
- Bed trolley
- Bar with trapeze
- Siderail pad
- Pump support
- V1100 mattress in 3 pieces

