




VF-300



VitalFlex ELITE LTC Bed



The VitalFlex Elite Homecare Bed is an all-in-one solution for users looking for **ergonomics, safety** and **style**. Its design-forward, quality construction is ideal for home and nursing environments. The VitalFlex Elite features Trendelenburg, Position, Cardiac Chair and a pedal break.



VitalFlex ELITE LTC Bed

Product Code: VF-300

Load capacity: 520 lbs / 235 kg

User Capacity: 440 lbs / 200 kg

Bed Weight: 233 lbs / 106 kg

Min.Height: 9.4 inches / 24 cm

Max Height: 31.5 inches / 80 cm

Max. Backrest Angle: 70°

Max. Legrest Angle: 30°

Sleep Deck Dimensions: 35.4" x 78.7" / 90 cm x 200 cm

Overall Dimensions: 41" x 84" / 104 cm x 214 cm

Mattress Guides: YES

Trendelenburg Position: BOTH WAYS

Trapeze: OPTIONAL



Folding Rail System

Half-rails can be quickly lowered and folded for safe patient access to the bed

Features & Benefits

- Tool-free installation
- Includes one full rail and two folding half rails
- Trendelenburg and Reverse Trendelenburg included
- Full wood cover panels
- Mattress guide at head and footrest
- 4-section steel-slat platform base
- Mattress guide at head and footrest
- All-metal sections are powder coated for durability
- 6-position mechanically adjusted footrest
- Four casters with central pedal break
- Pre-set Comfort Position
- 24V power supply for extra safety
- Ergonomic handset with locking key

Warranty: All parts covered for 24 months from date of purchase



Built to Last

Durable wood paneling and powder-coated steel are easy to maintain & ideal for long term care settings



Safe Positioning

Mattress guides keep sleep surface in place.



Comfort Position

Pre-set for comfortable TV viewing - ideal for cardiac and respiratory patients



Class-Leading Brakes

A convenient central foot brake system keeps bed in place with one simple motion