



# Jazzy<sup>®</sup>

## ELITE HD

The Jazzy<sup>®</sup> Elite HD offers heavy-duty construction, a 450 lbs. weight capacity, and excellent performance in a stylish, highly maneuverable package.



**Pride<sup>®</sup>**  
Mobility Products Co.  
— LIVE YOUR BEST<sup>®</sup> —

## SPECIFICATIONS

<b>Weight capacity</b>	450 lbs.*
<b>Maximum speed</b>	Up to 4.2 mph
<b>Ground clearance</b>	3"
<b>Turning radius</b>	24"
<b>Overall length</b>	42.5"
<b>Overall width</b>	24.5"
<b>Drive wheels</b>	14" solid
<b>Front wheels</b>	3" solid anti-tips
<b>Rear wheels</b>	6" solid casters
<b>Drivetrain</b>	Two-motor, front-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	60 amp controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 5 amp (standard)
<b>Per-charge range</b>	Up to 9.8 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	NF-22
<b>Battery weight</b>	37 lbs. each
<b>Base weight</b>	98.5 lbs.
<b>Standard seat weight</b>	51.5 lbs. (high-back)
<b>Max. seat dimensions</b>	24" x 24"
<b>Warranty</b>	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

\* Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18-in deep seat configurations.

## FEATURES

- Two-motor, front-wheel drive
- 14" knobby tires and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- 5 amp, off-board charger

## OPTIONS

- High-back seat with solid seat pan
- Synergy® seat
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder


## RECOMMENDED LIFTS

**Pride®**  
**Lifts & Ramps**

- Backpacker® Plus
- Outlander
- Backpacker® MV



## COLOR CHOICES



Red

Blue

**PRIDE®**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST® —

St. Catharines, Ontario • 888-570-1113 • [www.pridemobility.com](http://www.pridemobility.com)